

Scalable storage solutions to grow with your needs.

Infinity™ Medical Image Server



CODONICS

A new dimension in image storage

Infinity™ C/LX

Medical Image Server


The Infinity Medical Image Server combines simplicity and speed in a highly-reliable DICOM storage device. All Infinity models provide flexibility for permanent and temporary storage of DICOM images, reports and raw data in a compact unit. The unique, expandable design grows as your storage needs increase. Additionally, Infinity provides optional web viewing capabilities so medical images can be accessed over the network from any PC. Whether you are large institution looking for supplemental storage or a small facility transitioning to digital imaging, Infinity is the right choice.



Infinity C is the perfect low-cost storage solution for clinics and small hospitals as an alternative to PACS

Radiologist or Referring Physician

Infinity LX is ideal for facilities that need fast local storage to accommodate large data sets or as a complement to PACS




Features / Benefits

- Low-cost solution for storage of DICOM images and reports
- Designed to fill the gap between workstation storage and PACS
- Permanent and temporary image storage provides flexibility for multiple applications
- Remote network-based image viewing from a web browser*
- Compatible with existing DICOM viewers and PACS workstations
- Expandable as your storage needs increase**
- Removable media backup for offsite disaster recovery**
- Easy set up for non-technical users
- High performance architecture streamlines workflow with a powerful Oracle database and high-speed data transfers of up to 50MB/second (DICOM-store)
- Infinity models offer 3 TB to 8 TB of usable and scalable storage
- Compatible with all Codonics products including:
 - Virtua® Medical Disc Publisher for CD and DVD image distribution
 - Integrity® Medical Image Importer for reading and reconciling medical studies imported from disc
 - Horizon® Multi-media Dry Imager for hardcopy output to patients and referring physicians

*Optional feature
**Work-in-progress

Remote Image Viewing



Call Codonics today at **1-440-243-1198** or visit www.codonics.com for more information.

All registered and unregistered trademarks are the property of their respective owners. Specifications subject to change without notice. Printed in the U.S.A. Copyright © 2009 Codonics, Inc. 4/2009

17991 Englewood Drive
Middleburg Heights, OH 44130 USA
(440) 243-1198
(440) 243-1334 Fax
Email info@codonics.com
www.codonics.com