


eFilm Lite Viewer from Codonics

The eFilm Lite Viewer incorporates window and level presets, synchronized stacking, cine function, and advanced capabilities like volume rendering.



Codonics Virtua™ Medical Disc Publisher is the new standard for speed, efficiency and ease of use in automatic disc recorders. Virtua allows you to include the industry's leading viewers while its brilliant, full color disc labels increase your brand awareness.

For PACS and
general radiology

- Cinemode capable
- Window leveling, magnification, and measurements tools
- Annotations
- D.S.A. (digital subtraction angiography) to improve contrast
- Icon-based GUI
- MPR/MIP (3D) including stereo display mode

For more information contact Codonics at
800-444-1198 or **440-243-1198** and ask for
our Virtua Sales Department.

Contact your Codonics representative for viewer configurations and availability. Specifications are subject to change without notice. Contact your Codonics representative for the latest information. Copyright © 2007 Codonics, Inc.



17991 Englewood Drive
Middleburg Heights, OH 44130 USA
(440) 243-1198
(440) 243-1334 Fax
Email info@codonics.com
www.codonics.com